

KPL HistoMark[®] BLUE

For Localization of Alkaline Phosphatase-Labeled Reagents

Catalog No.
5510-0037 (55-70-00)

Size
500 mL

DESCRIPTION

KPL HistoMark[®] BLUE Substrate System is designed for visualization of alkaline phosphatase-labeled (AP) reagents. KPL HistoMark BLUE is a Fast Blue stain and KPL Contrast RED is a nuclear counterstain formulated from nuclear Fast Red. The substrate system provides a blue specific stain with red counterstain for immunohistochemical staining or immunoblotting.

KIT COMPONENTS

	Catalog No.	Volume
KPL PhThaloBLUE	5510-0039 (71-00-03)	10 mL
KPL Activator Solution	5570-0002 (71-00-01)	10 mL
KPL Buffered Substrate	5570-0003 (71-00-04)	50 mL
KPL Contrast RED	5540-0001 (71-00-05)	50 mL

Sufficient reagents are supplied to prepare 500 mL Substrate Solution (approximately 1000 slides).

STORAGE/STABILITY

- Reagents are stable for a minimum of one year stored at 2-8°C.
- Store KPL Contrast RED Solution tightly capped at room temperature.
- Discard KPL PhThaloBLUE Solution if black precipitate develops.
- Discard KPL Activator Solution or KPL Buffered Substrate Solution if yellow color develops.
- Warm all reagents to room temperature (24-28°C) before use.
- If a light precipitate is visible in KPL Buffered Substrate Solution, warm for 10 – 15 minutes in 37°C waterbath. Mix thoroughly by inversion until completely in solution.

REAGENTS NOT INCLUDED

1. Primary antibody.
2. AP-labeled reagents
3. Isopropyl alcohol.
4. Aqueous or xylene-based mounting media.
5. 0.1 M Tris-HCl (see **BUFFER PREPARATION**)
6. 1 M Citric Acid Free Acid (see **BUFFER PREPARATION**)

PREPARATION

- Substrate Solution (prepare immediately before use in Step 10)

NOTE: Prior to preparation, if a light precipitate is visible in KPI Buffered Substrate Solution, warm for 10 – 15 minutes in 37°C waterbath. Mix thoroughly by inversion until completely in solution.

- a. Add 0.5 mL KPL Buffered Substrate Solution to 5 mL reagent quality water.
- b. Mix 0.1 mL KPL PhThaloBLUE Solution with 0.1 mL KPL Activator Solution in a separate tube. Mix gently and allow to stand 3 minutes.
- c. After 3 minutes combine solutions from steps a and b. Mix thoroughly and use immediately.

- KPL Contrast RED Solution is supplied at use dilution.

PROCEDURE

1. Rehydrate paraffin embedded sections through graded alcohol (3 minutes each in 100%, 80%, 40% and 20% EtOH) to water. Other samples do not require rehydration. Frozen sections must be thoroughly dried before use.
2. KPL HistoMark BLUE reagents contain levamisole to block endogenous phosphatase activity. If additional blocking is required, apply Bouin's Solution or 1M citric acid 1 - 10 minutes.
3. Rinse slide 5 minutes in reagent quality water.
4. Soak in 0.1 M Tris-HCl for 3 - 10 minutes.

NOTE: Inorganic phosphate inhibits alkaline phosphatase activity. Avoid use of PBS or any solution containing phosphates.

5. Treat sample with primary antibody diluted in Tris-HCl for 15 - 20 minutes.

NOTE: Extended incubation may improve sensitivity.

6. Wash sample with Tris-HCl for 10 minutes.
7. Incubate sample with biotin-labeled antibody, directed against the primary antibody host species, 15 - 20 minutes. If using an AP-labeled secondary antibody, proceed to Step 9.
8. Wash as in Step 6.
9. Shake off excess buffer and incubate sample with AP-Streptavidin or AP-labeled secondary antibody diluted in Tris-HCl 15 - 20 minutes.

KPL HistoMark® BLUE

For Localization of Alkaline Phosphatase-Labeled Reagents

Catalog No.

5510-0037 (55-70-00)

Size

500 mL

10. Wash as in Step 6. (Prepare KPL Substrate Solution during this step.)
11. Shake off excess buffer and cover section with KPL Substrate Solution.
12. Incubate 10 minutes at room temperature out of direct light.
13. Rinse slide 2 - 3 minutes in reagent quality water.
14. Counterstain in KPL Contrast RED Solution 5 - 10 minutes.
15. Rinse thoroughly in reagent quality water until excess stain is removed from slide.
16. Air dry and mount in aqueous mounting medium. DO NOT USE XYLENE BASED MOUNTING MEDIA.

RESULTS

- Sites of enzyme activity range from pale to deep blue. Nuclei appear a contrasting pale pink to red.
- Sections not reacted with primary antibody as a negative control should not develop any blue tint.

NOTES

1. Always incorporate appropriate positive and negative controls.
2. Use substrate reagents immediately after mixing.
3. Instant development of blue color indicates that the primary antibody or phosphatase-labeled reagent must be further diluted.
4. Prolonged incubation in substrate may increase background and inhibit nuclear counterstaining.

BUFFER PREPARATION

0.1 M Tris-HCl

- a. Dissolve 121 g Tris in 500 mL reagent quality water.
- b. Adjust pH to 7.6 with 2 M HCl (approximately 300 mL).
- c. QS to 1 L with reagent quality water to obtain a 1 M stock.
- d. Dilute 1 part stock from Reagents Section, Step 5, with 9 parts reagent quality water and mix well.

1 M Citric Acid Free Acid

- a. Dissolve 192 g of citric acid free acid in 500 mL reagent quality water.
- b. QS to 1L with reagent quality water.

PRINCIPLE

The application of antibodies and other reagents such as avidin, streptavidin, etc., covalently coupled to calf intestine alkaline phosphatase in immunohistology is well documented (1, 2). The procedure described in this insert employs a simultaneous capture azo-dye technique, providing the research laboratory a method for precise localization of alkaline phosphatase-labeled reagents (3, 4). Primary aryl amines, when reacted with alkyl nitrites in acid media, form azo compounds (5). These react with substituted naphthols to produce highly chromogenic insoluble dyes. In this procedure the phosphate ester of 6-bromo-2-hydroxy-3-naphthoic acid (KPL Buffered Substrate Solution) is employed as the substrate. Enzymatic hydrolysis, in the presence of diazotized 4'-amino-2-,5-diethoxybenzanilide (KPL PhThaloBLUE Solution) results in the formation of a brilliant blue reaction product. Endogenous enzyme is eliminated by the incorporation of levamisole (6). It should be noted that a levamisole-resistant alkaline phosphatase has been demonstrated in some malignant cells from serous effusions (7). Additional blocking measures may be required (8, 9).

REFERENCES

1. Boorsma, D.M. (1984). *Histochemistry* 80: 103.
2. Jablonski, E., Moomaw, E.W., Tullis, R.H. et al. (1986). *Nucleic Acids Res.* 14: 6115.
3. Collings, L.A., Poulter, L.W., Janossy, G. (1984). *J. Immunol. Meth.* 75: 227.
4. Janckila, C.J., Yam, L.T., Li, C.Y. (1985). *Amer. J. Clin. Pathol.* 84: 476.
5. Burstone, M.S. (1962). *Enzyme Histochemistry and Its Application in the Study of Neoplasms*, Academic Press, New York, 88.
6. Van Belle, H. (1972). *Biochim. Biophys. Acta.* 289: 158.
7. Yam, L.T., Janckila, A.J., Epremian, B.E. et al. (1989). *Amer. J. Clin. Pathol.* 91: 31.
8. Pickel V.T., Joh, T.H., Reis, D.J. (1976). *J. Histochem. Cytochem.* 24: 792.
9. Molin, S.O., Nyrgen, H., Dolonius, L. (1978). *J. Histochem. Cytochem.* 26: 412.

The product listed herein is for research use only and is not intended for use in human or clinical diagnosis.